





# OPTIMISE COST AND QUALITY OF **GLAUCOMA** SURGERY

Our innovative surgical instrument solution provides total peace of mind, confidence and satisfaction.



**HASA OPTIX**  
Belgium

# Build your own sets for **GLAUCOMA** surgeries



## PREMIUM QUALITY

Premium quality disposable recyclable instruments offering constant quality with increased convenience and efficiency.

## COST CONTROL

HASA OPTIX instruments provide the ideal procurement solution by aiding the precise calculation of surgery costs.

## SIMPLIFY

Our pre-packed pouches help reduce the number of SKUs, contributing towards a simpler, more efficient repository solution.

## PEACE OF MIND


Our innovative solution provides total peace of mind ensuring the correct surgical instruments are easily accessible.

## INFECTION CONTROL

Our single use and pre-sterilized instruments ensure all surgeons can be confident with their patient safety of care.

## ENVIRONMENTALLY FRIENDLY

Supplied in recyclable pre-packed pouches relevant to each ophthalmic procedure, helping to reduce waste and environmental impact.




**HASARECYCLE**  
EXCELLENCE PROGRAM




At HASA OPTIX, we have decided to create a fully integrated re-use circle so that the stainless steel, our instruments are made of, can have a second life as non-medical objects.

## Anatomy Forceps




## Anatomy Speculums



All customized set and individual instrument are available in a box of 10, unless MOQ requirements are specified.

## Anatomy Scissors



**Sharp / Pointed**  
for tissue piercing and cutting



**Rounded**  
for cutting only.



**Straight**  
for when direct access is available or when a straight cut is required



**Curved**  
to reach more difficult areas or when a curved cut is required



**Angled**  
enables blades to be positioned at the correct angle to the issue, particularly useful when excising.



**Angled / Curved**




**Spring**  
enhances manual grip with excellent control of the tips for delicate procedures



**Bow or Ribbon**  
for non-delicate excisions such as extraocular (significant hand, wrist, or forearm movement may be required)

## Anatomy Needle Holders



**Fine / Micro**  
for finest sutures



**Standard**  
for general purpose



**Heavy**  
for larger sutures



**Straight\***




**Curved\***

\*Choice depends on users preference although curved shafts may give better visualisation of the wound.



**Flat**  
for enhanced manual grip and control of the tip (significant hand, wrist, or forearm movement may be required)



**Round**  
enhances micro control of the tips with minimum movement by rolling the handle between the fingers



**Lock Option**  
used to lock the jaws in the closed position, securing the needle during suturing

## HASA OPTIX's BESTSELLERS

### Forceps



551001




551187



551175



551199



551122



551131

- 551001 Toothed Forceps 0.12 Flat Handle
- 551187 Colibri Forceps 0.12mm Serrated Handle
- 551175 Tying Forceps Curved Jaws Flat Handle
- 551199 McPherson Lens IOL Forceps
- 551122 Notched Forceps Straight
- 551131 Colibri Notched Forceps Curved

### Scissors



653034



653055




553014



553016

- 653034 Westcott Scissors Sharp Curved
- 653045 Westcott Scissors Blunt Straight
- 653055 Westcott Scissors Blunt Curved
- 653056 Westcott Scissors Sharp Straight
- 553014 Vannas Curved Sharp Blades
- 553016 Vannas Straight Sharp Blades

### Needle Holders




552027



552022



552025



552005

- 552027 6mm Curved Jaws Round Knurled Handle
- 552022 10mm Curved Jaws Round Knurled Handle with Lock
- 552025 9mm Curved Jaws Round Knurled Handle
- 552005 6mm Straight Jaws Round Knurled Handle

### Marker



558003

- 558003 Scleral Marker - Double-ended 3.5/4.0mm

### Punch



554020

- 554020 Trabeculectomy Punch 0.75mm
- 554021 Trabeculectomy Punch 1.0mm
- 554022 Trabeculectomy Punch 1.5mm



Cataract

Cornea

Refractive

Glaucoma

Oculoplastic

Vitreoretinal

## SIMPLIFYING OPHTHALMIC SURGERY



Contact us to find your nearest distributor


@ info@hasaoptix.com

www.hasaoptix.com

+32.2.524.63.88

523 Avenue Louise  
1050, Brussels, Belgium

Follow us



**HASARECYCLE**  
EXCELLENCE PROGRAM

